


ECG Electrode Belt

Universal ECG Belt



Multi-use ECG belt that connects to all ECG devices

- 12-lead
- Compatible with All ECG Brands
- Fast & Easy to Place
- One Size Fits All - Elastic
- For Men & Women
- Reusable
- No Shaving, No Pain, No Disposables



- Connects to lead wires using Snap / Clip connectors
- 6 built in Ag / AgCl electrodes
- Provided with:
 - 4 limb clamps Ag / AgCl electrodes
 - Banana to Clip adaptors
 - Soft Carry Case



(for 3 & 4 mm lead wires)

Hygienic Protection ECG Belts




ECG Belts come with the integrated Sanitized® antimicrobial hygiene function for optimal protection.

What Does It Do?




The Sanitized® antimicrobial silver protects the ECG Belt. The silver additive is a scientifically undisputed substance, whose effect on bacteria is well known. The Sanitized® technology uses a glass ceramic capsule. This ensures that the antibacterial protective treatment is firmly anchored to the ECG Belt and thus achieves long-lasting protection.

Silver ions cannot diffuse from the ECG Belt thanks to the glass ceramic



The silver ion locates the bacteria via the moisture and deactivates it.



This is the effect silver has on bacteria:

1. Cell membrane is destabilized
2. Respiration is prohibited
3. Food (nutrient) intake is impeded
4. Cell division is inhibited

