

NORAV MEDICAL CONNECTIVITY

INTERFACES

		NEMS-Q	NEMS-A	NEMS-Web	PC-ECG Rest	PC-ECG Stress	Mobile ECG	NH-301 Holter	NECG	NBP-24NG ABPM	NBP-One ABPM	NSPIRO	NM-700 Telemetry
DICOM	Worklist				X	X ⁽¹⁾		X ⁽¹⁾					X ⁽¹⁾
	Report				X	X ⁽¹⁾		X ⁽¹⁾	X	X ⁽¹⁾	X ⁽¹⁾	X ⁽¹⁾	
HL7	TCP/IP Socket	X ⁽²⁾										X ⁽²⁾	
	Shared Folder	X	X		X ⁽¹⁾	X ⁽¹⁾		X ⁽¹⁾		X ⁽¹⁾	X ⁽¹⁾	X ⁽¹⁾	
	Web API				X ⁽⁴⁾								
GDT		X	X	X	X	X		X		X	X	X	X
FTP; SFTP ⁽¹⁾ ; FTPS ⁽⁴⁾				X ⁽⁴⁾	X ⁽¹⁾	X ⁽¹⁾	X ⁽⁴⁾	X ⁽¹⁾	X	X ⁽¹⁾	X ⁽¹⁾	X ⁽¹⁾	X ⁽¹⁾
E-mail					X ⁽³⁾	X ⁽³⁾	X	X ⁽³⁾			X		X ⁽³⁾
Amazon AWS S3 ⁽⁴⁾								X					
Web Upload ⁽⁴⁾					X	X		X					
Patient List (*.ini file)		X	X	X	X	X		X			X		X
NEMS-Q database		X			X	X	X	X	X	X	X	X	X
NEMS-A database			X		X	X		X	X	X	X	X	X
NEMS-Web database				X	X	X	X	X		X	X		X

DATA EXPORT

		NEMS-Q	NEMS-A	NEMS-Web	PC-ECG Rest	PC-ECG Stress	Mobile ECG	NH-301 Holter	NECG	NBP-24NG ABPM	NBP-One ABPM	NSPIRO	NM-700 Telemetry
DICOM	Encapsulated PDF				X	X ⁽¹⁾		X ⁽¹⁾		X ⁽¹⁾	X ⁽¹⁾	X ⁽¹⁾	
	Secondary Capture Image				X								
	12-Lead ECG Waveform				X				X				
HL7 ORU^R01 report		X			X	X ⁽¹⁾		X ⁽¹⁾		X ⁽¹⁾	X ⁽¹⁾	X ⁽¹⁾	
PDF report		X	X	X	X	X	X	X		X	X	X	X
XML format report		X	X			X	X			X		X	
FDA XML 12-Lead ECG					X				X				X
SCP 12-Lead ECG					X				X				
McKesson 12-Lead ECG					X								
MFER 12-Lead ECG					X		X						
RAW data ECG file					X	X	X	X ⁽⁵⁾					
Matlab format ECG data file					X								
Plain text report					X	X							
GDT report					X	X		X	X	X		X	X
Picture format report					X	X							
Excel file								X ⁽⁶⁾		X			

Note ⁽¹⁾ – by NoravScheduler;

Note ⁽²⁾ – optional, requires an additional fee;

Note ⁽³⁾ – POP3 e-mail client is required;

Note ⁽⁴⁾ – by Norav Uploader;

Note ⁽⁵⁾ – for DL800/DL900 records only;

Note ⁽⁶⁾ – for NR records only;