

NORAV MEDICAL CONNECTIVITY

INTERFACE		NEMS-Q	NEMS-A	NEMS-Web	PC-ECG Rest	PC-ECG Stress	Mobile ECG	NH-301 Holter	NECG	NBP-24NG ABPM	NBP-One ABPM	NSPIRO	NM-700 Telemetry
DICOM	Worklist	X ⁽²⁾		X ⁽²⁾	X	X		X ^(5,6)	X	X ^(5,6)	X ^(5,6)		X ⁽¹⁾
	Report	X ⁽²⁾		X ⁽²⁾	X	X		X ^(5,6)	X	X ^(5,6)	X ^(5,6)	X ⁽¹⁾	
HL7	TCP/IP Socket	X ⁽²⁾		X ⁽²⁾									
	Shared Folder	X ⁽²⁾		X ⁽²⁾	X ⁽¹⁾	X ⁽¹⁾		X ^(5,6)		X ^(5,6)	X ^(5,6)	X ^(5,6)	
	Web API				X ⁽⁴⁾								
GDT		X	X		X	X				X	X	X	X
FTP; SFTP ⁽¹⁾ ; FTPS ⁽⁴⁾					X ⁽¹⁾	X ⁽¹⁾	X		X	X ⁽¹⁾	X ⁽¹⁾	X ⁽¹⁾	X ⁽¹⁾
E-mail					X ⁽³⁾	X ⁽³⁾	X	X ⁽³⁾					X ⁽³⁾
Amazon AWS S3 ⁽⁴⁾					X ⁽⁴⁾	X ⁽⁴⁾		X ⁽⁴⁾					
Web Upload ⁽⁴⁾					X ⁽⁴⁾	X ⁽⁴⁾		X ⁽⁴⁾					
Patient List (*.ini file)		X	X		X	X							X
NEMS-Q database		X			X	X	X	X	X	X	X	X	X
NEMS-A database			X		X	X		X	X	X	X	X	X
NEMS-Web database				X	X	X	X	X	X	X	X		X

DATA EXPORT		NEMS-Q	NEMS-A	NEMS-Web	PC-ECG Rest	PC-ECG Stress	Mobile ECG	NH-301 Holter	NECG	NBP-24NG ABPM	NBP-One ABPM	NSPIRO	NM-700 Telemetry
DICOM	Encapsulated PDF	X		X	X	X		X ^(5,6)	X	X ^(5,6)	X ^(5,6)	X ⁽¹⁾	
	Secondary Capture Image	X		X	X	X							
	12-Lead ECG Waveform	X		X	X				X				
HL7 ORU^R01 report		X ⁽²⁾		X ⁽²⁾	X	X ⁽¹⁾		X ^(5,6)		X ⁽¹⁾	X ⁽¹⁾	X ⁽¹⁾	
PDF report		X	X	X	X	X	X	X	X ⁽²⁾	X	X	X	X
XML format report						X	X					X	
FDA XML 12-Lead ECG					X				X				X
SCP 12-Lead ECG					X				X ⁽²⁾				
McKesson 12-Lead ECG					X								
MFER 12-Lead ECG					X		X						
RAW data ECG file					X	X	X						
Matlab format ECG data file					X								
Plain text report					X	X							
GDT report					X	X			X	X		X	X
Picture format report					X	X							
Excel file					X	X							
EDF format ECG data file								X					

Note ⁽¹⁾ – by NoravScheduler;

Note ⁽²⁾ – optional, requires an additional fee;

Note ⁽³⁾ – POP3 e-mail client is required;

Note ⁽⁴⁾ – by Norav Uploader;

Note ⁽⁵⁾ – by NEMS-Q;

Note ⁽⁶⁾ – by NEMS-Web;