





Norav Medical

Welcome to Norav Medical, a global leader in ECG solutions and products, delivering exceptional quality and innovation since the nineties. For over 30 years, Norav Medical has set the standard for ECG excellence with our state-of-the-art technologies and products.



Norav Medical is a leading medical device manufacturer specializing in the development and production of advanced cardiovascular monitoring systems. Since the nineties, we have provided innovative and reliable solutions for healthcare professionals and patients worldwide. Our products are designed to improve patient outcomes, reduce healthcare costs, and enhance the overall quality of care. We are committed to excellence in every aspect of our business, from product development to customer service.

Thank you for choosing Norav Medical, where our legacy of innovation and excellence in ECG solutions continues to lead the industry.

Our offerings are distinguished by their diagnostic accuracy, comprehensive features, and user-friendly interfaces. We emphasize reliable connectivity to EHR/HIS systems, including HL7 and DICOM standards. Many Norav products offer wireless connectivity to smartphones, tablets, and PCs, supporting efficient workflows for healthcare professionals. We serve large healthcare facilities, clinics, hospitals, OEMs, and medical research establishments. Our expert software and hardware customization teams provide fast-response support, ensuring our solutions meet the unique needs of our diverse clientele.

NORAV FOUR LAYER SOLUTION

1

Test Devices

- All types of ECG tests
- High-precision ECG Devices
- Portable and Stationary Models
- Reliable Signal Acquisition

2

Diagnostic Applications

- Advanced Algorithms
- User-friendly
 Interfaces
- Real-time Analysis

3

Seamless Connectivity

- Integration with EHR/HIS Systems
- Wireless Connectivity (Smartphones, Tablets, PCs)
- HL7 and DICOM Standards

4

ECG Management

- Unified Platform for Data
 Management (Web/local)
- View, Analyze & Manage ECGs
- Streamlined Workflows
- Improved Patient Care











RELIABILITY

Designed for robust operation in diverse and challenging environments, with both wired and wireless acquisition.



FLEXIBILITY

In-house software development enables prompt response and custom solutions for market needs.



ACCURACY

Providing the best ECG signal quality and automatic measurement and interpretation.



QUALITY

Operating under the most rigorous standards, we are committed to serving the medical industry.



UNIVERSAL CONNECTIVITY

Solutions for HL7, GDT, and DICOM connectivity, as well as custom HIS interfaces.



ECG SPECIALIZATION

With over 30 years of experience in developing, manufacturing, and innovating ECG PC-based and ECG systems.



DURABILITY

Long-lasting products for heavy use, guaranteed by highquality manufacturing technology and excellent raw materials.



SAFETY & EFFICIENCY

Maintaining safety and essential performance of the electrocardiograph system under safety and EMC regulations.

OUR PORTFOLIO



Resting ECG



Stress ECG



Holter ECG



ABPM






Spirometry



Resting ECG systems are fundamental tools in cardiovascular diagnostics, designed to capture the electrical activity of the heart when a patient is at rest. Norav Medical's Resting ECG solutions are renowned for their accuracy, reliability, and user-friendly interfaces. These devices are essential for detecting various cardiac conditions, from arrhythmias to myocardial infarctions, providing healthcare professionals with crucial insights for patient care.

NORAV'S RESTING ECG SOLUTION









Resting ECG Key Features



High Diagnostic Accuracy

Advanced signal processing ensures precise measurements even during physical exertion.



User-Friendly Design

Intuitive interfaces and automated features simplify operation.



Comprehensive Connectivity

Integration with EHR/HIS systems (HL7, DICOM) ensures seamless data management.



Robust and Durable

Built to withstand rigorous use in clinical and fitness settings.



Enhanced Workflow Efficiency

Streamlined data management and reporting capabilities improve clinical workflows.



Mobile Application

Real-time ECG monitoring and data sharing on-the-go for enhanced accessibility.



Diagnostic Software

Powerful tools for in-depth analysis and interpretation of ECG data.



ECG Management System

A unified platform to display, analyze, and manage all types of ECG tests, providing a comprehensive solution for healthcare facilities.

NEMS-WEB: ECG MANAGEMENT SYSTEM

Web-based Interface, Multi Sites & Users

Serial ECG Comparison


Unlimited Archive and History

ECG Acquisition and Reporting

Real-Time Resting ECG Access

Advanced Resting ECG Analysis

NEMS-Web seamlessly integrates with Resting ECG systems, providing real-time data access and advanced analysis for precise







ECG Mobile Application Resting ECG

ECG's Anytime, Anywhere

Transform your mobile device into a powerful diagnostic electrocardiograph with the Norav Mobile ECG App. Designed specifically for Resting ECG tests, this app enhances accessibility and provides healthcare professionals with real-time ECG monitoring and analysis capabilities, no matter where they are.





Mobile ECG Acquisition

Perform and record resting ECGs directly on your Android, iOS, or Windows device, turning your smartphone or tablet into a fully functional ECG machine.



Real-Time ECG Transmission

Instantly transmit ECG results as PDF attachments for a second opinion or directly to the emergency room, ensuring quick and informed medical responses.



Data Comparison on-the-Go

Easily store and access previous ECG tests, enabling side-by-side comparisons of historical data for more accurate diagnostics.



Efficient and Portable

Whether in the clinic or remote locations, Norav Mobile ECG offers a reliable way to monitor patients, ideal for mobile healthcare, home visits, and emergency care.



Built-in MEANS Rest ECG Interpretation

Integrated with an advanced automatic diagnosis system, offering reliable interpretations and enhanced diagnostic confidence.



12-Lead ECG Recording

Capture high-quality, diagnosticgrade 12-lead ECGs directly on your mobile device, ensuring precise readings.





iOS and Android Compatibility

The app works seamlessly across both iOS and Android platforms, enabling widespread accessibility on various devices.



Derived Leads

Offers 5-lead ECGs that can be derived into 12-lead readings on Android and Windows devices for flexible use (optional feature).



Comprehensive Display and Printout

Offers automatic display and printout options during tests, making data sharing and documentation efficient.



1200M Resting ECG

Classic & Affordable ECG

The 1200M is an affordable yet highly effective ECG device, balancing cost and performance. Ideal for general practice, community clinics, and low-resource settings.



12

International Standards FDA. CE



EHR / HIS Interface HL7, DICOM, GDT



Connection

Wired USB



Powering

5V (USB)



Our benefits:

- Cost-effective
- High-Resolution ECG
- Accuracy at High Sampling Rate
- MEANS interpretation
- Straightforward Installation and User Friendly
- On-Line Filtering & Automatic Baseline Correction
- Automatic interpretation
- Multiple languages

- Caliper QRS and ECG Measurements
- Comprehensive display and printout
- Physician remarks and signature
- Email ECG results directly from your Mobile, Tablet or PC
- Automatic print during test
- Network abilities and data sharing
- Laser or ink printout
- Bulk export ECG files to XML, Excel, TXT

High-Resolution ECG

The 1200HR provides high-resolution ECG recordings, essential for detailed cardiac analysis. Suitable for cardiology clinics, research facilities, and advanced diagnostic centers.

1200HR Resting ECG



12 , 16



AnalysisLate Potentials,



International Standards FDA, CE



EHR / HIS Interface HL7, DICOM, GDT,



Connection

Wired USB



Powering

5V (USB)

- 24-Bit Resolution and High Sampling Rate
- Signal Averaging (Late Potentials) and Heart Rate Variability applications
- Clean & Stable ECG Traces
- MEANS interpretation
- Straightforward Installation and User Friendly
- On-Line Filtering & Automatic Baseline Correction
- Automatic interpretation
- Advanced analysis tools

- Caliper QRS and ECG Measurements
- Comprehensive display and printout
- Physician remarks and signature
- Email ECG results directly from your Mobile, Tablet or PC
- Automatic print during test
- Network abilities and data sharing
- Multiple languages
- Laser or ink printout
- Bulk export ECG files to XML, Excel, TXT







NR-1207-3 Resting ECG

All-in-One ECG

The NR-1207-3 combines multiple functionalities in a single device. offering comprehensive ECG solutions. Suitable for hospitals, specialized cardiology units, and emergency departments.



International **Standards** 3,12 FDA, CE



EHR / HIS Interface HL7, DICOM



Connection

Wireless Bluetooth



Powering

AA batt



Our benefits:

- Rest, Stress, Holter, and Telemetry
- Small, Lightweight, and Water Resistant
- **MEANS** interpretation
- 5 electrode cable derived to 12 Leads in Android & MS Windows (Optional)
- Accurate ECG Recordings
- High-Resolution ECG Processing
- MS Windows is Compatible with all ECGs

- iOS, Android, and MS Windows Compatible for Rest ECG
- **Automatic interpretation**
- Caliper QRS and ECG Measurements
- · Comprehensive display and printout
- Physician remarks and signature
- Email ECG results directly from your Mobile, Tablet or PC
- Multiple languages
- Laser or ink printout
- Bulk export ECG files to XML, Excel, TXT

Mobile Wireless ECG

The NR-1207-E is designed for wireless operation, providing mobility and convenience for ECG diagnostics. Ideal for home healthcare, mobile clinics, and emergency response teams.

NR-1207-E Resting ECG



3,12



International Standards FDA, CE



EHR / HIS Interface HL7, DICOM



Connection

Wireless Bluetooth



Powering

AA batt

- Small, Lightweight, and Water Resistant
- MEANS interpretation
- 5 electrode cable derived to 12 Leads in Android & MS Windows (Optional)
- Accurate ECG Recordings
- High-Resolution ECG Processing
- Advanced analysis tools
- MS Windows is Compatible with all ECGs

- iOS, Android, and MS Windows Compatible for Rest ECG
- Automatic interpretation
- Caliper QRS and ECG Measurements
- Comprehensive display and printout
- · Physician remarks and signature
- Email ECG results directly from your Mobile, Tablet or PC
- Multiple languages
- Laser or ink printout
- Bulk export ECG files to XML, Excel, TXT





NECG SE 1200 Pro

Advanced 12-Lead Electrocardiograph

The SE-1200 Pro provides high-quality 12-lead ECG recordings with Glasgow automatic interpretation, including DICOM connectivity, designed for clinics and hospitals.







EHR / HIS Interface DICOM, HL7, GDT





Our benefits:

- 80 kHz High Sampling Rate for Precision
- Glasgow automatic interpretation
- Auto Pacemaker Detection with Independent Detection Channel
- Gesture-Based Amplification and Measurement of Waveforms
- Advanced Al Algorithms for Disease Detection and Notifications
- Seamless Integration with EHR/HIS (DICOM, HL7, GDT)
- Automatic Signal Recording with Reverse Electrode and Lead Connection Warnings
- Serious Disease Red Alerts for Urgent Diagnoses

- Comprehensive Report Comparison Tools
- Averaged Template Analysis with Manual Recalculation
- Pacemaker Detection
- Self-Adapt AC Filtering for Clear ECG Signals
- 8-Inch High-Resolution with Color Display
- Touch-Screen Interface
- Dedicated Function Keys
- Compact and Lightweight Design
- Built-In Thermal Printer
- Durable, Easy-to-Clean
- Battery Backup
- Multi-Language Support



Conventional ECG

The NECG SE-301 is a state-of-the-art 3-lead ECG device, offering unparalleled accuracy and ease of use. Suitable for use in emergency departments, clinics, and general practice settings.

NECG SE 301 Resting ECG





800 ECGs

International Standards FDA



EHR / HIS Interface HL7, DICOM



Connection
Wireless
WiFi



Batter\



*Memory can be extended using a USB drive or memory card up to 8 GB

- Full touch operation
- Clean & Stable ECG Traces
- Seamless Integration with HIS/HER
- Up to 5-hours battery life
- Internal memory up to 800 ECGs *NEMS is not included
- Easy transfer of stored data to NEMS
- Straightforward Installation and User Friendly

- Anti-noise technology with comprehensive filters and high CMRR
- Automatic ECG measurements and interpretation
- High-quality digital filters for correction of baseline shift, signal noise from EMG, and high frequency
- Continuous Impedance Test with Alarm Display







ECG BELT Resting ECG

Time Savings, Universal. No shave, No gel, No Disposable Electrodes

The ECG Belt offers a novel approach to ECG monitoring, ensuring quick and accurate electrode placement. Suitable for outpatient clinics, general practice, and research studies.







Size Adjustable , Pediatric







*Pediatric belt available for children and smaller patients

Our benefits:

- Designed for 12-lead Rest ECG
- Easy application with adjustable fit
- Secure electrode placement for accurate results
- Comfortable for long-term use •
- Reduces setup time, improves comfort
- Reusable, cost-efficient design
- Fits all body sizes
- Compatible with all 12-lead ECG brands

- Durable and long-lasting materials with Hygienic Protection
- Includes soft carry case with protection
- built-in 6 precordial electrodes
- No shave, no gel needed
- Snap or clip lead wire connectors
- Ideal for clinics, practices, and research
- Consistent electrode positioning for accuracy
- Quick and easy to place

NECG Rolling Carts

The NECGTRolling carts with adjustable height and storage, perfect for the NECG Conventional Electrocardiographs, provides a convenient and organized solution for housing ECG equipment.


NECG CARTResting ECG



- Compatible with all NECG line,
 NECG SE-1200 Pro and NECG SE-301
- Caster base plate with 5 wheels
- Adjustable height
- Basket attachment
- Optional cable arm attachment available








Stress ECG systems are essential for assessing cardiac function under physical stress. These systems help detect ischemia, arrhythmias, and other cardiac abnormalities that may not be present at rest. Norav Medical's Stress ECG solutions are designed for accuracy, reliability, and ease of use, providing healthcare professionals with critical information to evaluate and manage cardiac health.

NORAV'S STRESS ECG SOLUTION









Stress ECG Key Features



High-Resolution ECG

Ensures accurate ST-measurements and arrhythmia detection with clean traces.



Advanced Signal Quality

Provides clear, stable ECG signals during stress tests.



Seamless System Integration

Compatible with HIS/EHR and CPET systems for easy data management.



Robust and Durable

Built to withstand rigorous use in clinical and fitness settings.



Real-Time Data Playback

Allows reviewing missed episodes post-test for accurate diagnosis.



Wireless Patient Mobility

Enables patient movement without signal loss during testing.



Comprehensive Diagnostic Software

Features Duke and HRR Scores for in-depth analysis.



Customizable Test Protocols

: Supports standard and userprogrammed stress test options.

NEMS-WEB: ECG MANAGEMENT SYSTEM

Web-based Interface, Multi Sites & Users

Enhanced Security & Data Protection


Unlimited Archive

ECG Acquisition and Reporting

Seamless Integration with HIS/EHR via HL7 & DICOM

Patient History and Worklist

NEMS-Web seamlessly integrates with Stress ECG systems, providing real-time data access and advanced analysis for precise







1200S Stress ECG

Classic Stress ECG

The 1200S is a classic stress ECG system known for its reliability and ease of use in stress testing environments. Ideal for general practice, community clinics, and hospitals.



Lead Systems

12



Score

HRR & Duke Treadmill



International Standards FDA, CE



EHR / HIS Interface HL7, DICOM, GDT



Connection

Wired USB



Powering

5V (USB)



Our benefits:

- High-Resolution ECG Processing
- Clean & Stable ECG Traces
- Seamless Integration with HIS/EHR
- Straightforward Installation and User-Friendly
- Standard/Unlimited Custommade, as well as Treadmill Test (TMT) Stress Test Protocols
- Automatic Arrhythmia Detection
- Automatic arrhythmia detection, print, and capture
- Duke Treadmill Score

- HRR Heart Rate Recovery Score
- ST/HR Index
- Clean Trace Synthesized ECG Rhythm Preserving Actual Heart Rate, Rhythm, and Beat Shapes
- Full disclosure enables postprocessing of saved data
- Interfaces with a variety of CPET systems
- Automatic print during test
- Plain or thermal paper printing options
- Automatic control of a variety of treadmills and ergometer

High-Resolution ECG - Echo Ready

The 1200HR-T delivers high-resolution ECG recordings with superior accuracy, featuring high sampling rate and Echo Ready Trigger Out functionality. Ideal for advanced cardiac diagnostics in research and clinical settings



Lead Systems

12



Bit Resolution

24



International Standards FDA, CE

EHR / HIS Interface [HL7, DICOM, GDT]

1200HRT Resting ECG



Connection

Wired USB



Powering

5V (USB)

- Lead-Off Detection
- Pacemaker Detection
- Echo Ready Trigger OutClean & Stable ECG Traces
- Playback of Episodes
- Standard/Unlimited Custommade Protocols for Ergometers, as well as Treadmill Test (TMT) Stress Test Protocols
- Seamless Integration with HIS/EHR
- Straightforward Installation and User-Friendly
- Selectable Sync Lead
- Late Potentials & Heart Rate
 Variability Application Compatible

- Automatic arrhythmia detection, print, and capture
- Duke Treadmill Score
- HRR Heart Rate Recovery Score
- ST/HR Index
- Full disclosure enables postprocessing of saved data
- Interfaces with a variety of CPET systems
- Automatic print during test
- Plain or thermal paper printing options
- Automatic control of a variety of treadmills and ergometers
- Selectable Sync Lead* *With selected models





1200W

Stress ECG

Ultimate Wireless Stress ECG

The 1200W is a top-of-the-line RF wireless stress ECG system that delivers high-performance diagnostics in a wireless format. Ideal for cardiology clinics, stress testing labs, and hospitals.



12





International Standards FDA, CE



EHR / HIS Interface HL7, DICOM, GDT



Connection

Wireless RF



Powering

AA batt



Our benefits:

- Accurate ECG Recordings
- High resolution ECG processing
- Clean & stable ECG traces
- Playback of episodes
- Echo-ready wireless acquisition
- Standard/unlimited custommade, as well as Treadmill Test (TMT) Stress Test protocols
- Seamless integration with HIS/EHR
- Straightforward installation and user-friendly
- Selectable sync lead

- Automatic arrhythmia detection, print, and capture
- Duke Treadmill Score
- HRR Heart Rate Recovery Score
- ST/HR Index
- Standard as well as unlimited user-programmed Stress Test protocols
- Clean Trace Synthesized ECG Rhythm Preserving Actual Heart Rate, Rhythm, and Beat Shapes
- Interfaces with a variety of CPET systems
- Automatic control of a variety of treadmills and ergometers

Blood Pressure Monitor

The Tango M2 is a specialized blood pressure monitor designed for use during stress testing. Ideal for cardiology clinics, stress testing labs, and research facilities.

Our benefits:

- Automated Blood Pressure Measurements
- For cardiac Stress and Exercise Testing
- 7" color LCD display
- Accurate and reliable blood pressure measurements during stress tests
- Designed specifically for treadmill and ergometer stress testing
- Automatic BP measurements, even during patient movement
- Artifact filtering for clean, accurate readings

- Compatible with all major stress test systems
- BP readings unaffected by noise or patient motion
- User-friendly interface for easy operation
- Compact and durable design for clinical use
- Quick and efficient measurements, minimizing interruptions during stress testing
- Seamless integration with ECG and stress systems
- USB Port for Easy Upgrade

Tango M2 Stress ECG



USB, RS232

Powering

AC mains





NR-1207-3 Stress ECG

All-in-One ECG

The NR-1207-3 combines multiple functionalities in a single device, offering comprehensive ECG solutions for stress testing. Suitable for hospitals, specialized cardiology units, and stress testing labs.



International
Standards
FDA, CE



EHR / HIS Interface HL7, DICOM



Connection
Wireless
Bluetooth




AA batt

3,12

NR NORAV





- Small, Lightweight, and Water Resistant
- Accurate ECG Recordings
- High-Resolution ECG Processing
- Playback of episodes
- Automatic arrhythmia detection, print, and capture
- Duke Treadmill Score
- long-lasting products with robust designs

- HRR Heart Rate Recovery Score
- ST/HR Index
- Standard as well as unlimited user-programmed Stress Test protocols
- Clean Trace Synthesized ECG Rhythm Preserving Actual Heart Rate, Rhythm, and Beat Shapes
- Interfaces with a variety of CPET systems
- Automatic control of a variety of treadmills and ergometers
- Automatic print during test

Mobile Wireless ECG

The NR-1207-E is designed for wireless operation, providing mobility and convenience for stress ECG diagnostics. Ideal for home healthcare, mobile clinics, and stress testing labs.

NR-1207-E Stress ECG







EHR / HIS Interface HL7. DICOM



Connection

Wireless Bluetooth



AA batt

- Rest, Stress
- Small, Lightweight, and Water Resistant
- Accurate ECG Recordings
- High-Resolution ECG Processing
- Playback of episodes
- Automatic arrhythmia detection, print, and capture
- Duke Treadmill Score
- long-lasting products with robust designs

- HRR Heart Rate Recovery Score
- ST/HR Index
- Standard as well as unlimited user-programmed Stress Test protocols
- Clean Trace Synthesized ECG Rhythm Preserving Actual Heart Rate, Rhythm, and Beat Shapes
- Interfaces with a variety of CPET systems
- Automatic control of a variety of treadmills and ergometers
- Automatic print during test







Treadmills Stress ECG

Medical-Grade Treadmills

Medical-grade with robust design for heavy duty treadmills designed for accurate and reliable exercise testing. Ideal for stress testing labs, cardiology clinics, and rehabilitation centers.







Our benefits:

- Heavy-duty design with triple-thickness for maximum durability and comfort
- Medical-grade treadmills specifically designed for cardiac stress testing
- Customizable protocols tailored to patient-specific needs
- Easy-to-use programmable controller for stress test management (selected models)
- Nonslip rubber coating for enhanced patient safety and stability
- Versatile, supporting a range of stress test requirements
- Seamless integration and control with Norav's Stress Test ECG devices and software
- High energy absorption for smooth, comfortable operation during tests
- Precision-engineered for accurate stress test results
- Available in a variety of models to suit different clinical needs
- High energy absorption for smooth, comfortable operation during tests

Medical-Grade Ergometers

Medical-grade with robust design for heavy duty ergometers designed for accurate and reliable exercise testing. Ideal for stress testing labs, cardiology clinics, and rehabilitation centers.

Our benefits:

- Heavy-duty reclined ergometers for robust performance
- Precise load control throughout the entire load range
- Universal interfaces for seamless communication with external devices
- All key data displayed on a single control terminal for easy monitoring
- ECG connection for accurate measurement and recording during tests
- Adjustable seat for enhanced patient comfort and support
- Remote control for convenient adjustment of test positions and seat settings
- Available in a variety of models to meet different clinical needs
- Fully controlled by Norav's Stress Test software

Ergometers Stress ECG












Holter ECG systems are essential tools for continuous cardiac monitoring, providing valuable data over extended periods to detect intermittent arrhythmias and other cardiac anomalies that might not be captured during a standard ECG test. Norav Medical's Holter ECG solutions are designed for high diagnostic accuracy, patient comfort, and ease of use, making them ideal for both clinical and outpatient settings.

NORAV'S HOLTER ECG SOLUTION









Holter ECG Key Features



High Diagnostic Accuracy

Advanced algorithms ensure precise detection of cardiac events.



Continuous Monitoring

Provides detailed recordings over extended periods (up to 14 days).



Comprehensive Connectivity

Seamless integration with EHR/HIS systems (HL7, DICOM).



Compact and Lightweight

Designed for patient comfort during long-term monitoring.



Respiration Trends

Monitors patient breathing patterns during the recording.



Audible Diary Recording

Captures patient-reported events with voice memos.



Patient Physical Activity Monitoring

Tracks the patient's physical activity alongside ECG data.



Customizable Reporting

Tailored reports for detailed analysis of long-term data.

NEMS-WEB: ECG MANAGEMENT SYSTEM

Web-based Interface, Multi Sites & Users

Enhanced Security & Data Protection


Patient Hookup and Study Upload

Report Review and Confirming

Seamless Integration with HIS/EHR via HL7 & DICOM

Patient History and Worklist

NEMS-Web seamlessly integrates with Holter ECG systems, providing real-time data access and advanced analysis for precise







Holter ECG systems are essential tools for continuous cardiac monitoring, providing valuable data over extended periods to detect intermittent arrhythmias and other cardiac anomalies that might not be captured during a standard ECG test. Norav Medical's Holter ECG solutions are designed for high diagnostic accuracy, patient comfort, and ease of use, making them ideal for both clinical and outpatient settings.

Features	NR-302	NR-314	NR-1207	NR-1207-3	NR-314-P
3 Channels ECG Recording	/	/	/	/	/
12-lead ECG Recording			/	/	
Pacemaker detection	/	✓	/	✓	✓
Respiration (TIV) recording		/	/	/	
Accelerometry sensor		/	/	/	/
Voice note recording		/	/	✓	
USB data transfer		/	/	/	/
Bluetooth Wireless Transmission		/	✓	✓	✓
Holter recording * wireless ECG at the same time				✓	✓
Recording time (up to)	48 H	336 H	336 H	336 H	336 H

NR-301

The ultimate Holter Analysis solution

- Intuitive Best-in-Class Editing Tools
- 24 hr 336 hr Recording
- QT variability
- Heart Rate Variability
- HL7 Interface
- Remote Analysis / Uploader









NR-302 Holter ECG

Affordable Holter

The NR-302 is an affordable Holter ECG device designed to balance cost and performance for continuous cardiac monitoring. Ideal for general practice, community clinics, and low-resource settings.



REC •

Recording

Up to 48h



International Standards FDA, CE



EHR / HIS **Interface** HL7, DICOM, GDT



Connection

Wired USB



Powering

AA batt



Our benefits:

- Small, Lightweight, and Water Resistant
- 250 Hz ECG sampling rate
- Arrhythmia, QT, ST, HRV, and Pacemaker analysis
- Atrial Fibrillation detection algorithm
- Automatic interpretation and report generation
- Physical activity trends tracking.
- Lead-off and pacemaker detection

- Flashcard reader and USB for data transfer
- 512MB memory for extended recordings
- On-line filtering and baseline correction
- Display and automatic print during tests
- Physician remarks and signature support
- Remote analysis and upload options
- Water-resistant with IP22 rating
- Voice recording for patient-reported events

14-Days Holter

The NR-314 provides up to 14 days of continuous monitoring, offering detailed and long-term cardiac data. Suitable for outpatient clinics, cardiology practices, and research studies.

NR-314 Holter ECG



REC •

Recording

Up to 14 Days

International Standards

FDA, CE

EHR / HIS **Interface** HL7, DICOM, GDT



Connection USB. Bluetooth



AA batt

- Small, Lightweight, and Water Resistant
- Supports 3-lead, 5-lead, and 12-lead
- High sampling rates: 250, 500, 1000 Hz
- Arrhythmia, QT, ST, HRV & Pacemaker analysis
- Atrial Fibrillation detection algorithm
- Automatic interpretation and report generation
- Lead-off detection for accurate monitoring
- Pacemaker detection for improved diagnostics

- Respiration monitoring via thorax impedance
- Tracks physical activity trends
- 2GB memory for extended storage
- Flashcard and USB data transfer
- Voice recording for patient events
- Online filtering and baseline correction
- Display and automatic print during tests
- Supports physician remarks and signature
- Bluetooth ECG transmission during patient hookup
- Water-resistant (IP22) and durable



NR-1207 Holter ECG

Ultimate Holter

The NR-1207 is a versatile Holter ECG device offering comprehensive monitoring and advanced analysis features. Ideal for cardiology clinics, hospitals, and research facilities.



REC •

Recording

Up to 14 Days



International Standards FDA, CE



EHR / HIS Interface HL7, DICOM, GDT



Connection USB, Bluetooth



Powering

AA batt

Our benefits:



- Supports 3-lead, 5-lead, and 12-lead
 ECG configurations
- High sampling rates: 250, 500, 1000 Hz
- Arrhythmia, QT, ST, HRV, & Pacemaker analysis
- Atrial Fibrillation detection algorithm
- Automatic interpretation and report generation
- Lead-off and pacemaker detection
- Respiration monitoring via thorax impedance
- Physical activity trends tracking

- 2GB memory for extended storage
- Flashcard and USB data transfer
- Voice recording for patient events
- Online filtering and baseline correction
- Bluetooth ECG transmission during patient hookup
- Water-resistant (IP22) and durable
- Supports physician remarks and signatures
- Comprehensive display and automatic print during tests

All-in-One ECG

The NR-1207-3 combines multiple functionalities in a single device, offering comprehensive ECG solutions for Holter monitoring. Suitable for hospitals, specialized cardiology units, and research facilities.

NR-1207-3 Holter ECG



3,12

REC •

Recording

Up to 14 Days

International Standards FDA, CE



EHR / HIS Interface HL7, DICOM, GDT



Connection
USB,
Bluetooth



Powering

AA batt

- All-in-one ECG solution: Rest, Stress, Holter, Telemetry
- Compact design for patient comfort
- Supports 3-lead, 5-lead, and 12-lead ECG configurations
- High sampling rates: 250, 500, 1000 Hz
- Arrhythmia, QT, ST, HRV & Pacemaker analysis
- Atrial Fibrillation detection algorithm
- Automatic interpretation and report generation
- Lead-off and pacemaker detection

- Respiration monitoring via thorax impedance
- · Physical activity trends tracking
- 2GB memory for extended storage
- Flashcard and USB data transfer
- Voice recording for patient events
- Online filtering and baseline correction
- Bluetooth ECG transmission fduring patient hookup
- Water-resistant (IP22) and durable
- Supports physician remarks and signature
- Comprehensive display and automatic print during tests





NR-314-P Holter ECG

Compact size 14-day Holter The NR314-P is a compact and lightweight Holter ECG device

designed for up to 14 days of continuous monitoring. : Suitable for outpatient clinics, cardiology practices, and research studies.



REC • Recording Up to 14 Days



International EHR / HIS **Standards Interface** FDA, CE HL7, DICOM, GDT





Connection USB. Bluetooth



Powering Internal Lithium-ion polymer battery













Our benefits:

- Lightweight: only 41g for enhanced patient comfort
- Long battery life with built-in Lithium-polymer battery
- Supports 3-lead, 5-lead, and derived 12-lead ECG configurations
- High sampling rates: 128, 256, 512, 1024 Hz
- Arrhythmia, QT, ST, HRV & Pacemaker analysis
- Atrial Fibrillation detection algorithm
- Automatic interpretation and report generation
- Lead-off and pacemaker detection
- Tracks physical activity trends
- 2GB memory for extended storage
- Flashcard and USB data transfer
- Bluetooth ECG transmission
- Water-resistant (IP64) and durable
- Supports physician remarks and signatures
- Comprehensive display and automatic print during tests

Fits Seamlessly Into Your patient's Life

Up to 14 days of continuous ECG recording allows you to go about your regular activities uninterrupted, . The lightweight and discreet design makes the device virtually invisible.

Precision Where It Matters

The Norav NR-314-P adheres directly to the body using medical-grade gel electrodes, not a patch. This ensures the highest quality signal acquisition without interference for incredibly accurate and reliable test results.

Peace of Mind Through Smart Analysis

The Holter software* with an advanced Atrial Fibrillation algorithm provides precise AFib detection and analysis. Combined with arrhythmia, QT, ST segment, and heart rate variability assessments, you get a comprehensive analysis of cardiac health and function.

Designed for Efficiency

Quick integration with HIS/EHR systems along with DICOM compatibility streamlines workflows. Self-customised printouts, email sharing of results, and multi-language options simplify use.





Norav Medical's Telemetry Rehab systems are essential for continuous cardiac monitoring during rehabilitation, offering real-time data to ensure patient safety and enhance recovery protocols. With arrhythmia, ST, and SpO2 monitoring, alarms, and comprehensive reporting, our solutions provide unmatched diagnostic accuracy. Complete 12-lead surveillance and clean exercise ECG traces ensure precise monitoring even in challenging rehabilitation settings, making Norav's systems the ideal choice for clinical and home rehabilitation, designed for ease of use and patient comfort.

NORAV'S TELEMETRY REHAB SOLUTION



All in One ECG

Telemetry ECG System

Exercise SpO2


Medical-Grade Ergometers


Medical-Grade Treadmills

Norav Medical's Rehabilitation ECG System offers complete 12 lead telemetry surveillance, enabling real-time display and storage of traces.

With PC-based central station, the system offers comprehensive reports of all patient activity, while keeping operation simple.

Security
Compliance
Diagnostics
Data storage











Telemetry Rehab Key Features



High Diagnostic Accuracy

Advanced algorithms ensure precise detection of cardiac events.



Continuous Monitoring

Real-time cardiac data ensures patient safety during rehabilitation.



Comprehensive Connectivity

Seamless integration with EHR/HIS systems (HL7, DICOM).



Patient Comfort

Lightweight, portable designs for ease during physical rehabilitation.



Treadmill & Ergometer Control

Direct control of exercise devices during sessions.



SpO2 Monitoring

Optional feature for tracking oxygen saturation



Simultaneous Test Control


Manage and monitor multiple rehabilitation tests at the same time for enhanced efficiency.



Advanced ECG Management

A unified platform to display, analyze, and manage rehabilitation ECG data.

REAL-TIME DETECTION, ALARMS AND REPORTING.







NR-1207-3 Telemetry

All-in-One ECG

The NR-1207-3 combines multiple functionalities in a single device, offering comprehensive ECG solutions for telemetry and rehabilitation. Suitable for hospitals, specialized cardiology units, and rehabilitation centers.



Lead Systems
3,12
International
Standards
FDA, CE



EHR / HIS Interface HL7, DICOM



Connection

Wireless Bluetooth



Powering

AA batt



Our benefits:

- Rest, Stress, Holter, and Telemetry •
- Arrhythmia, ST, and SpO2 monitoring, alarms and reporting.
- Small, lightweight, and waterresistant design (IP22)
- Clean exercise ECG traces
- High sampling rates: 250, 500, 1000 Hz
- Seamless integration with HIS/EHR systems
- Six Minute Walk Test (6MWT) and
 Shuttle Walk Test support
- Treadmill and Ergometer control

- Real-time wireless ECG transmission via Bluetooth
- High-resolution ECG processing
- Physical activity trends and patient monitoring
- Supports monitoring of multiple patients simultaneously
- Simple, intuitive operation and easy installation
- Online filtering and baseline correction for clean signals
- 2GB memory for full disclosure backup

Telemetry Rehab ECG

The NR-314-T is a compact telemetry ECG device designed for continuous monitoring during cardiac rehabilitation. Suitable for outpatient clinics, rehabilitation centers, and home healthcare.

NR-314-T Telemetry







EHR / HIS Interface HL7, DICOM



Connection
Wireless
Bluetooth



Powering

AA batt

- Small, lightweight, and waterresistant design (IP22)
- Seamless integration with HIS/EHR systems
- Arrhythmia, ST, and SpO2 monitoring, alarms and reporting.
- Clean exercise ECG traces
- Supports Six Minute Walk Test (6MWT) and Shuttle Walk Test
- Simple operation with intuitive user interface
- Control of treadmills and ergometers

- Wireless transmission via Bluetooth for remote analysis
- Advanced report generation
- Online filtering and automatic baseline correction
- High-resolution ECG processing for accurate results
- Tracks physical activity trends during rehab
- patient-reported events
- Tracks physical activity trends during rehab





Treadmills[®] Telemetry

Medical-Grade Treadmills

Medical-grade treadmills designed for accurate and reliable exercise testing and rehabilitation. Ideal for rehabilitation centers, cardiology clinics, and stress testing labs.









Our benefits:

- Heavy-duty design with triple-thickness for maximum durability and comfort
- Medical-grade treadmills specifically designed for cardiac stress testing
- Customizable protocols tailored to patient-specific needs
- Easy-to-use programmable controller for stress test management (selected models)
- Nonslip rubber coating for enhanced patient safety and stability
- Versatile, supporting a range of stress test requirements
- Seamless integration and control with Norav's Stress Test ECG devices and software
- High energy absorption for smooth, comfortable operation during tests
- Precision-engineered for accurate stress test results
- Available in a variety of models to suit different clinical needs
- High energy absorption for smooth, comfortable operation during tests

Medical-Grade Ergometers

Medical-grade ergometers designed for accurate and reliable exercise testing and rehabilitation. Ideal for rehabilitation centers, cardiology clinics, and stress testing labs.

Our benefits:

- Heavy-duty reclined ergometers for robust performance
- Precise load control throughout the entire load range
- Universal interfaces for seamless communication with external devices
- All key data displayed on a single control terminal for easy monitoring
- ECG connection for accurate measurement and recording during tests
- Adjustable seat for enhanced patient comfort and support
- · Remote control for convenient adjustment of test positions and seat settings
- · Available in a variety of models to meet different clinical needs
- Fully controlled by Norav's Stress Test software

Ergometers elemetry









NSPO2 Telemetry

Exercise SPO2 with BLE

The nSpO2 02 is an exercise SpO2 monitor designed to measure blood oxygen levels during physical activity. Ideal for rehabilitation centers, cardiology clinics, and stress testing labs.







Battery LifeUp tp 44 hours



Our benefits:



- Fast and accurate SpO2 and pulse rate readings
- · Continuous monitoring during stress tests and rehabilitation
- PureSAT® SpO2 technology for precise accuracy
- · Compact, wrist-worn design for patient comfort
- Large display with enhanced sensor connector
- Long-lasting battery life (up to 44 hours)
- Pulse-by-pulse averaging algorithm for precision
- 1,080 hours of memory at a 4-second sampling rate
- · Ideal for small, portable oximetry during exercise
- SpO2 range: 0% to 100%, pulse rate range: 18-321 BPM
- Wireless data transmission for real-time monitoring
- Durable and reliable for high-intensity activities

Telemetry Monitoring System

The NM-700 is a state-of-the-art monitoring system designed for continuous telemetry during cardiac rehabilitation. Ideal for cardiac rehabilitation centers, hospitals, and home healthcare.


NM-700Telemetry

- Complete 12 leads surveillance
- Continuous telemetry monitoring for cardiac and pulmonary rehab
- Clean exercise ECG traces
- High-resolution ECG processing for clean, accurate traces
- Real-time Arrhythmia, ST, and SpO2 detection, alarms and reporting.
- Supports 3, 6, and 12-lead ECG configurations
- PC-based central station for real-time monitoring
- Advanced algorithms for accurate arrhythmia detection
- Treadmill and ergometer control
- Online filtering and automatic baseline correction
- Seamless integration with HIS/EHR and PACS/DICOM
- · Remote analysis and uploader for easy data sharing
- Simple operation with user-friendly installation
- PDF reports for efficient data review
- Comprehensive display and automatic print during tests
- Long battery life and easy-to-use design















Norav TeleEVENT is a state-of-the-art Mobile Cardiac Telemetry (MCT) solution designed for healthcare providers who demand precision, reliability, and improved patient outcomes. Featuring the NR-314-P wearable chest strap, a user-friendly mobile app, and the secure NEMS-Web management platform, Norav TeleEVENT streamlines continuous ECG monitoring, rapid event detection, and real-time reporting—empowering clinicians to diagnose faster and intervene earlier.

NORAV'S TELEEVENT SOLUTION

From Patient to Provider—A Simple, Three-Step Process

Norav TeleEVENT enables patients to remain at home while their heart is monitored continuously and discreetly. Physicians gain access to immediate, actionable data, reducing diagnostic delays and enhancing care efficiency.

Connect



ECG Device + Chest Belt

- High-Precision ECG Recording
- Lightweight, comfortable NR-314-P chest strap
- · Minimal patient intervention required
- Up to 10 days continuous monitoring battery
- Up to 30 days record of full disclosure
- Recording of up to 500 events
- Holter and dual-mode capabilities
- An allergy-free, sticker-free monitoring belt


Transmit



Patient App

- Automatic and manual wireless ECG data transmission,
- Automatic identification and transmission of critical arrhythmias.
- Secure HIPAA-compliant communication

Review & Respond



Web-Based ECG Management System

- Customizable Notifications configurations
- Centralized patient data management
- Comprehensive reporting







WHY HEALTHCARE PROVIDERS CHOOSE NORAV TELECARDIO?

Scalable & Flexible

Edeal for clinics, hospitals, and dedicated monitoring centers—adaptable to various clinical environments.



High-Precision ECG Analysis

Advanced algorithms ensure accurate detection of arrhythmias (pauses, AFib, bradycardia, tachycardia).

Fully Wireless Integration

Streamlined workflows with simple smartphone connectivity, minimizing patient burden.



medical

ECG Excellence

Real-Time Notifications & Reports

Instant notifications for critical events, supporting timely interventions.



Patient-Centric Design

Non-invasive, comfortable, and convenient—improving compliance and diagnostic yield.

Comprehensive Web-Based ECG Management


NEMS-Web supports meaningful use by receiving, managing, archiving, and exporting diagnostic test results via a web-based interface that is customizable, intuitive and user-friendly.

THE NORAV TELEEVENT MOBILE APP SERVES AS THE VITAL LINK BETWEEN PATIENT AND PROVIDER.
BEYOND TRANSMITTING ECG DATA, THE APP ENSURES A SMOOTH, USER-FRIENDLY EXPERIENCE THAT
ENHANCES COMPLIANCE AND SUPPORTS CONTINUOUS MONITORING WITH:

MOBILE APP FEATURES







Telemetry Rehab systems are crucial for continuous cardiac monitoring during rehabilitation, providing real-time data to ensure patient safety and optimize recovery protocols. Norav Medical's Telemetry Rehab solutions are designed for high diagnostic accuracy, ease of use, and patient comfort, making them ideal for clinical and home rehabilitation settings.

NORAV'S ABPM SOLUTION


Norav Medical's Ambulatory Blood Pressure Monitoring is a state-of-the-art motion tolerant monitor that offers accurate 24-hour ambulatory blood pressure monitoring.

By recording up to 250 / 300 measurements of systolic and diastolic blood pressure and heart rate, it highly assists clinicians with better diagnosis and

management of hypertension.













ABPM Key Features



High Diagnostic Accuracy

Advanced algorithms ensure precise blood pressure measurements.



Continuous Monitoring

Provides comprehensive blood pressure data over 24 hours or more.



Comprehensive Connectivity

Seamless integration with EHR/HIS systems (HL7, DICOM).



Interpretive Summary

Delivers clear, easy-to-understand summaries of monitoring results for clinical use.



Smart Inflator

Prevents over-inflation by considering previous measurements, enhancing patient comfort.



Motion Tolerance

Ensures accurate readings, even when the patient is in motion.



Standards Compliance

Ensures accurate readings, even when the patient is in motion.



Customizable Reporting

Tailored reports for detailed analysis of long-term data.

NEMS-WEB: ECG MANAGEMENT SYSTEM

Web-based Interface, Multi Sites & Users

Enhanced Security & Data Protection


Patient Hookup and Study Upload

Report Review and Confirming

Seamless Integration with HIS/EHR via HL7 & DICOM

Patient History and Worklist

NEMS-Web seamlessly integrates with ABPM systems, providing real-time data access and advanced analysis for precise











Holter Blood Pressure

The NBP One is a compact and reliable Holter blood pressure monitor designed for continuous 24-hour monitoring. Suitable for outpatient clinics, cardiology practices, and hypertension management programs.

REC •

Recording

24h

Orbit



International Standards

FDA, CE



EHR / HIS Interface HL7, DICOM, GDT

Connection


Wired USB



Powering

2x AA batt





- Stores up to 250 measurements for extended data collection
- Lightweight and compact for patient comfort
- Adjustable ABPM Orbit cuffs available in two sizes (2 and 3)
- BHS & ESH compliance for accuracy and safety
- Day/night sleep time recording for personalized insights
- Customizable protocols for hypertension management
- Precise control in hypertensive situations
- Flash memory for secure data storage
- Seamless integration with GDT, HL7, and HIS/EHR
- USB 3.0 (USB-C) connectivity for data transfer
- Clear, concise, and personalized patient reports
- Long battery life with AA batteries (over 250 measurements)
- Portable and easy to operate for ambulatory monitoring
- Seamless integration with HIS/EHR and ECG Management System (NEMS)

Holter Blood Pressure

NBP One

NORAV

The Mobil-O-Graph is a state-of-the-art Holter blood pressure monitor providing detailed and accurate data for comprehensive blood pressure assessment. Ideal for general practice, cardiology clinics, and hypertension management programs.

Mobil-O-Graph II Ambulatory BP

REC •

Recording

24h



International **Standards**

FDA. CE



EHR / HIS Interface HL7, DICOM, GDT



Powering

2x AA batt

Our benefits:

Accurate 24-hour monitoring of blood pressure and heart rate

Records up to 300 measurements for detailed analysis

Lightweight and compact design for patient comfort

ABPM cuff available in multiple sizes (S, M, L, and optional XS, XL)

Day/night sleep time recording for personalized insights

Bluetooth and USB connectivity for fast data transfer

Compliant with BHS and ESH requirements for hypertension management

Custom-made and default monitoring protocols

Precise control in hypersensitive situations

Seamless integration with HIS/EHR and ECG Management System (NEMS)

Simple installation and user-friendly operation

- Personalized patient reports for efficient diagnosis
- Long battery life, supporting over 300 measurements
- Clear and concise data collection and retrieval











Spirometry is a key diagnostic tool for assessing lung function and diagnosing respiratory conditions such as asthma, chronic obstructive pulmonary disease (COPD), and other pulmonary disorders. Norav Medical's Spirometer solutions are designed for high diagnostic accuracy, ease of use, and patient comfort, making them ideal for clinical and outpatient settings.

NORAV'S SPIROMETER SOLUTION

Norav's Spirometer is a fast and easy system that provides accurate and reliable Spirometry measurements.

Optimal for pulmonary analysis, the Spirometer System is the perfect tool for asthma management, early detection of acute respiratory disorders and for the use in pharmaceutical trials.














Norav Medical ECG Management System (NEMS-Web) is an inclusive data management solution that syndicates integrated studies and patient records in a single application.

In an era where precision and efficiency in patient care are paramount, Norav Medical's NEMS-Web emerges as a beacon of innovation in ECG solutions. Designed to seamlessly integrate with diverse healthcare systems across the globe, NEMS-Web offers unparalleled diagnostic accuracy, user-friendly interfaces, and Seamless connectivity options. Customized to meet the precise demands of healthcare professionals globally, Our platform ensures that every heartbeat tells a story, fostering better patient outcomes through technology. Discover the future of ECG monitoring with NEMS-Web, where every element is crafted with excellence and every connection opens a new horizon in healthcare.

View, Analyze, and Manage ECGs


the simple way

NEMS-Web supports meaningful use by receiving, managing, archiving, and exporting diagnostic test results via a web-based interface. Using NEMS-Web, you can quickly create and confirm all types of ECG tests, as well as spirometry and blood pressure tests. Through a standard



web browser, users can access with safe and secure connection from everywhere around the globe, to perform, edit, and view test results.











NEMS-Web Key Features



HL7 / DICOM Support

Advanced interoperability with universal medical data protocols



EHR / HIS Connectivity

Easily integrate diverse healthcare systems for enhanced global connectivity

Telehealth Access

from patients home or

the system using

Perform ECG examinations

remote clinics and access

Smartphone, Tablet or PC



Empower Patients

View/analyze multiple ECG

• Integrate patient-generated ECG

Trusted Clinical Diagnostics Automated accurate interpretation. Tools to enhance clinical decisions.

formats.

• Remote monitoring and access options

Actionable Insights

- Dashboards with timely notifications
- Tailored customization settings
- Custom real-time reporting Serial comparison tools

Flexible and Scalable

- Hardware and virtual options
- Migrate data easily across systems.
- Grow from clinics to large hospitals

Cost-Effective Healthcare

Delivery: Improved diagnostic processes and streamlined operations lead to significant cost savings for healthcare institutions.



Multivendor Support

Seamless data exchange across multi-device ecosystems, and third-party ECG equipment support.

Cloud/Premise Support

Access the server from the

or on Premise (VM)

Cloud (AWS/Azure/Google)



Robust Security

Comprehensive SSL, twofactor authentication for data protection, and digital signature.



Workflow Simplified

User-friendly user interfaces enhance workflow and patient care for optimized efficiency.



Effortless Installation & **Maintenance**

Easy server and client's distribution for windows platform

NEMS-Web System Highlights

Web-based Interface, Multi Sites & Users

Serial ECG Comparison

Unlimited Archive and History


ECG Acquisition and Reporting

Enhanced Security & Data Protection

REAL-TIME Alerts and Notifications


Notifications system

Third-party PDF reports import





NEMS-Web Cloud Scanning Center











About Us

With over 30 years of expertise in ECG systems.
Norav Medical excels in developing, manufacturing, and innovating PC-based solutions. We prioritize safety and efficiency, adhering to stringent safety and EMC regulations. Our global network ensures expert and responsive local service. We also provide customizable solutions to meet the unique needs of healthcare professionals.



Mission and Vision

Our mission is to enhance cardiac care through innovative ECG solutions that provide healthcare professionals with the tools they need for accurate diagnosis and effective treatment. We are committed to quality, reliability, and continuous improvement, ensuring our products meet and exceed industry standards.



Global Presence

Norav Medical, with regional headquarters in Delray Beach, Florida (USA) and Mainz-Kastel (Germany), markets its products through a global dealer network. The company is ISO13485-certified and FDA-listed, and its products are FDA- and CE-certified, and ensure compliance with safety and EMC regulations to maintain the safety and essential performance of our electrocardiograph systems.



Support and Services

Our expert software and hardware customization teams provide fast-response support to OEMs and clinical research establishments. We are dedicated to delivering excellent customer service and technical support to ensure our products meet the evolving needs of healthcare providers worldwide.

BECOME A NORAV DISTRIBUTOR

Join the Norav Medical Family

Are you ready to be part of a global leader in ECG solutions? As a Norav Medical distributor, you will have the opportunity to offer state-of-the-art ECG products and services renowned for their quality, accuracy, and innovation. We provide comprehensive support, training, and resources to ensure your success in promoting and selling Norav Medical products.

Take the next step in your business growth and become a part of the Norav Medical family. Contact us today to learn more about this exciting opportunity!
























CONTACT US

info@norav.com

VISIT OUR WEBSITE

www.noravmedical.com

