

NORAV HOLTER NR-314-P



The New Holter NR-314-P

14 Days of Continuous ECG in Just 41g!



Rechargeable
Lithium-ion
polymer




Lightweight
41g



Water Resistant
IP64

Experience the future of cardiac monitoring with the Norav Holter NR-314-P. This compact and lightweight device is designed to be virtually invisible, allowing for discreet and comfortable wear. Despite its small size, it packs a powerful recording with up to 14 days of continuous recording, ensuring comprehensive cardiac monitoring. With water splashing protection (IP64), the Norav Holter NR-314-P is built to withstand daily activities, providing reliable and accurate ECG readings no matter what your day brings. Choose the Norav Holter NR-314-P for a seamless, efficient, and unobtrusive cardiac monitoring experience.



2-electrodes
Single channel ECG



CONTACT US
info@norav.com



VISIT OUR WEBSITE
www.noravmedical.com

UNINTERRUPTED

Daily activities

MONITORING THAT FITS YOUR LIFE



Fits Seamlessly Into Your patient's Life

Up to 14 days of continuous ECG recording allows you to go about your regular activities uninterrupted. The lightweight and discreet design makes the device virtually invisible.



Precision Where It Matters

The Norav NR-314-P adheres directly to the body using medical-grade gel electrodes, not a patch. This ensures the highest quality signal acquisition without interference for incredibly accurate and reliable test results.



Peace of Mind Through Smart Analysis

The Holter software* with an advanced Atrial Fibrillation algorithm provides precise AFib detection and analysis. Combined with arrhythmia, QT, ST segment, and heart rate variability assessments, you get a comprehensive analysis of cardiac health and function.



Designed for Efficiency

Quick integration with HIS/EHR systems along with DICOM compatibility streamlines workflows. Self-customised printouts, email sharing of results, and multi-language options simplify use.



The Results You Can Trust

Water-resistant construction and proven accuracy exceed the highest industry standards. FDA and CE provide further validation you're getting the best quality monitoring.




Showering is allowed while wearing the device, thanks to its IP64 water resistance.



Avoid direct water spray on the sensor. It is recommended to shower with your back towards the showerhead



Do not submerge in water! Refrain from swimming or bathing while using the device.



*Supported by Holter analysis software NH-301 v 4.0 and up


A COMPREHENSIVE SOLUTION FOR HOLTER ECG TESTS

Norav Medical's Holter Analysis System is a solution that provides complete Holter ECG monitoring for up to 14 days. With an array of small versatile Holter recorders to choose from the NH-301 is the ultimate Holter Analysis solution. Retrieve, View and Print ECGs.


Acquire



Analysis



Storage & Share



NH-301

Ultimate Holter Analysis System

Step into a new era of cardiac monitoring with Norav's NH-301 software, a comprehensive solution designed for extended ECG monitoring of up to 14 days. The NH-301 software seamlessly integrates with HIS/EHR and PACS/DICOM, facilitating remote analysis and PDF report generation. Its high-resolution recording and advanced analysis capabilities precisely detect and analyze arrhythmias, QT, ST, HRV, and atrial fibrillation, providing a robust foundation for accurate diagnostics.

EVENT MANAGEMENT:

Sort events for enhanced diagnostic precision, now with more examples for Min/Max HR & RR.

DATA ANALYSIS:

Relabel beats with unique colors; artifact beats exclusion for increased reliability.

USER INTERFACE:

Zoomable trends, physical activity display, and dual view for selected beats enhance your analysis experience.

MORE CAPABILITIES:

Direct Excel export for tabular reports and continuous data backup




5-electrodes
3 channel Monitoring



Belt+electrode
3 channel



4-electrodes
3 channel ECG



5-electrodes
with extenders
3 channel monitoring



3-electrodes
3 channel ECG



Norav Medical

Norav Medical is a leading medical device manufacturer that specializes in the development and production of advanced cardiovascular monitoring systems. With over 30 years of experience in the medical industry, Norav Medical is dedicated to providing innovative and reliable solutions for healthcare professionals and patients worldwide. Our products are designed to improve patient outcomes, reduce healthcare costs, and enhance the overall quality of care. At Norav Medical, we are committed to excellence in every aspect of our business, from product development to customer service.



Resting ECG



Stress ECG



Holter ECG



ABPM



Telemetry Rehab



Spirometry

USA

601 N. Congress Ave. Unit 105, Delray
Beach FL 33445
Tel: (+1) 561-274-4242

GERMANY

Christof-Ruthof-Weg 10, D55252
Mainz-Kastel
Tel: (+49) 6134.567983-0

NORAV
medical
ECG Excellence