


SMALLER. SMARTER. STRONGER. COMING SOON...



NORAV'S NEWSLETTER

With the milestone of MDR CE certification now achieved, Norav Medical is moving confidently into a new era of growth and innovation. Our focus on compliance and quality has paved the way for the launch of the next generation 1200W v3 Stress ECG – lighter, more compact, and more powerful, designed to set a new benchmark in performance and usability.

This achievement marks only the beginning: with breakthrough technologies and upcoming product innovations already in development, we continue to expand the boundaries of what ECG diagnostics can deliver. Together with our global distributor network, we are ready to shape the future of cardiac care.

Product is important, but having great product distribution is more important.



LET'S MEET AT MEDICA 2025

MEDICA



At the upcoming MEDICA show, we will continue to be located in the same familiar spot as in previous years, but this time we will present an innovative and newly designed 50-square-meter booth. It's a great opportunity to meet and discuss future collaborations. You can schedule a meeting with your area's Sales Manager.

November 17-20, 2025
Hall 9, Booth #E33 Düsseldorf, Germany.



Hospitalar 2025

May 20-23, 2025, São Paulo

WHX Miami 2025

June 11-13, 2025, Miami

We had a remarkable presence at WHX Miami 2025 and Hospitalar 2025, two of the largest healthcare exhibitions in the Americas. Miami brought together over 16,000 healthcare professionals and more than 1,300 exhibitors from 110+ countries, while Hospitalar in São Paulo hosted 1,200 exhibitors and 70,000 visitors from 30+ countries. Both events offered an exceptional stage to showcase Norav's advanced ECG portfolio, including the NR Series and TeleEVENT solution under development. Most importantly, we were delighted to reconnect with our dedicated long-term distributors and to welcome new partners into the Norav family. These connections strengthen our network and set the stage for continued growth across North and Latin America.

Top image: the Norav team at the WHX Miami 2025 exhibition in Miami.

Bottom image: The Norav team at the HIMSS 2025 exhibition in Las Vegas.

HIMSS 2025

March 3-6, 2025, Las Vegas

HIMSS®

25

AT HIMSS 2025 IN LAS VEGAS, A PREMIER DIGITAL HEALTH EVENT ATTRACTING TENS OF THOUSANDS OF GLOBAL HEALTH IT LEADERS AND INNOVATORS, NORAV HIGHLIGHTED ITS COMMITMENT TO CONNECTIVITY AND REMOTE CARDIAC CARE WITH NEMS-WEB AND TELEEVENT. THE EXHIBITION PROVIDED AN INVALUABLE OPPORTUNITY TO EXCHANGE INSIGHTS WITH HEALTHCARE EXECUTIVES, IT LEADERS, AND CLINICIANS. WE WERE ESPECIALLY PLEASED TO MEET OUR LOYAL DISTRIBUTORS ONCE AGAIN AND TO GREET NEW PARTNERS WHO JOINED US AT THIS IMPORTANT EVENT.



**COMING
SOON**



ENHANCED 1200W STRESS ECG – LIGHTER, SMARTER, AND STRONGER

Get ready for the next generation of Norav's 1200W Wireless Stress ECG System – a new evolution in **performance, connectivity, and design**. The upcoming model introduces a more compact structure, stronger RF communication, and enhanced pacemaker detection, all while being lighter and easier to use than ever before. Engineered to deliver precise cardiac diagnostics with improved stability and usability, this system represents Norav's vision for the future of stress testing. Currently under final development, **THE NEW 1200W IS NOT YET AVAILABLE - BUT A FIRST LIVE PREVIEW WILL TAKE PLACE AT MEDICA 2025 IN DÜSSELDORF. WE INVITE YOU TO VISIT OUR BOOTH AND BE AMONG THE FIRST TO EXPERIENCE THE FUTURE OF WIRELESS CARDIAC DIAGNOSTICS.**

KNOW MORE:

- **Refined Wireless Performance**
New RF architecture delivers higher transmission power (8 dB) and stronger connection stability.
- **Lighter, Compact Design**
Reduced device weight (252 g) and a slimmer profile improve handling and patient comfort.
- **Advanced Pacemaker Detection**
Enhanced digital detection at 128 ksps ensures precise rhythm analysis for complex cardiac signals.
- **Simplified Power Management**
Operates with a single AA battery, providing extended usage with lower maintenance.
- **Optimized Internal Architecture**
A new single-board structure replaces the previous three-PCB design, improving reliability and efficiency.
- **Full Compatibility**
Works seamlessly with Norav's stress ECG ecosystem – ensuring smooth transition and consistent performance.



**AVAILABLE
NOW**



User manuals now in 6 languages

DISCOVER AND DOWNLOAD NORAV'S **NR SERIES** USER MANUALS IN 5 NEW LANGUAGES FROM THE DISTRIBUTOR PORTAL TODAY. IT'S AVAILABLE NOW AND READY TO SUPPORT YOUR BUSINESS.

DOWNLOAD



STAY UPDATED

FOLLOW US



in





**COMING
SOON**

MINISPIRO - THE NEW SPIROMETER: EXPANDING THE POWER OF PRECISION DIAGNOSTICS

We are excited to introduce the new Spirometer, an advanced respiratory testing solution that expands Norav Medical's diagnostic ecosystem beyond cardiac care.

Compact, accurate, and effortless to use, the Minispiro provides real-time pulmonary function testing directly on your PC through the Spiro Software. Designed to support healthcare professionals in hospitals and clinics, it enables fast, reliable measurements for both routine screening and advanced diagnostic applications.

Fully compatible with Norav's ecosystem and EHR/HIS connectivity standards, the Minispiro represents the next step in our mission to deliver comprehensive, integrated, and data-driven diagnostic solutions.

[!\[\]\(d66ff64371a51729ac8c1cdaa685ba6f_img.jpg\) **READ MORE**](#)

KNOW MORE:

- **Real-Time Testing**
Perform a full range of pulmonary tests (FVC, VC, MVV, pre/post bronchodilator) with real-time data visualization.
- **Automatic BTPS Conversion**
Built-in temperature sensor ensures automatic calibration and uncompromised accuracy.
- **Seamless EMR/EHR Connectivity**
Integration via HL7, GDT, FHIR, and Exchange protocols for smooth data sharing and patient record updates.
- **Pediatric Mode**
Interactive on-screen animations enhance collaboration and accuracy for pediatric patients.
- **Plug & Play Simplicity**
USB-powered and software-driven for instant use-no setup, no internal memory, no calibration required.
- **Hygienic & Flexible Operation**
Supports both reusable and disposable turbines, offering safety, efficiency, and low maintenance.
- **Certified & Compliant**
Conforms to international standards including MDR 2017/745, FDA, CE, and ATS/ERS 2019 guidelines.





HOLTER NH-301 (V.4.0.64)

- **Personalized Event Names & Enhanced Tagging** - Edit event number manually , rename and organize events for a more intuitive review experience.
- **Noise Exclusion Option** – remove noisy X-beats to ensure cleaner, more accurate HR/RR analysis.
- **Simplified Navigation & Uniform Printing** – move through events with Next/Previous buttons and print exactly as displayed.
- **Page Screen-** On the Event page, highlight the focused event. Add a filter for Events (Checkbox list of all events with the option to check/uncheck all).

SOFTWARE
UPDATES



STRESS ECG – V5987

- **New ST/HR Index** – an advanced diagnostic parameter for improved ischemia evaluation.
- **Live Playback from Last Event** – review ECG data instantly during testing.
- **Expanded Device Support** – now compatible with Power Bike ergometers.
- **Password-Protected Access** – new secure login for PC-ECG Stress and Rest systems



ABPM MODULE (NEMS)

- **Automated Interpretive Reports** – clear results aligned with JNC7 & AHA guidelines.
- **Dipper / Non-Dipper Status** – identifies hypertension patterns and provides actionable insights.
- **Enhanced Clarity** – sleep/awake averages and category color coding for faster clinical interpretation.



NEMS A, NEMS-Q (V2.7.8) & NEMS-WEB

- **NR-314-P Support** – full integration with “Check ECG” hookup mode.
- **DICOM Archive & Worklist** – enables seamless connectivity with NECG SE 1200 PRO.
- **Localized Interface** – added German language support for ABPM.
- **Improved Summary Control** – toggle automatic interpretations on or off per preference




NEMS Web!

UPDATES!

V1.0.1.0

**AVAILABLE
NOW**



ENHANCED WORKFLOW & REPORTING

- **Custom Summary Templates** – personalize Holter and ABPM report layouts with your logo and contact info.
- **Imported Records Dashboard** – monitor and manage incoming data with live status tracking.
- **Dual Holter Report Mode** – automatically generate both full and 24-hour analysis (enables dual billing).
- **Integrated PDF Viewer** – keyboard shortcuts and auto-fit viewing for faster analysis.



SYSTEM INTEGRATION & CONNECTIVITY

- **Extended HL7 ORU Output Holter** – includes detailed SVE/VES ectopy data, tachycardia episodes, and HR metrics.
- **Expanded DICOM Support** – ensuring seamless connectivity with Norav Holter and Stress systems.
- **Improved Security Policies** – stronger password enforcement with configurable length and complexity.



REAL-TIME NOTIFICATIONS & ALERTS

- **In-App, Email, and SMS Notifications** – configurable alert options for faster response.
- **Custom Recipient Groups** – send notifications to specific users or departments.
- **Continuous Patient Monitoring** – receive instant alerts for new or abnormal ECG results.



EFFICIENCY & ADMIN TOOLS

- **Flexible User Access** – administrator rights no longer required to run the client.
- **Maintenance from Browser** – perform system maintenance directly within NEMS-Web.
- **Improved Alert Management** – immediate warnings when patient data is accessed by multiple users.



2025 Price List

Make sure you have the most updated price list

MiniSpiro

- **MINI-SPIRO-USB-01**
Mini,USB Spirometer,Incl2 Flow turbines,carboard,blue pouch
- **SPIRO-FLOW-TURBINE-01**
Flow turbine,MIR,Box 60 pcs,for Minspir USB



NEW

NSpiro



OUT OF STOCK

The NSpiro device is now discontinued and unavailable for ordering. We offer the new Spirometer - MiniSpiro. Discover more about these advanced options at the [link](#), ensuring quality and satisfaction for your customers.

Wall Holders

- **HOLDER WALL-NR-01**
Holder,Wall,NR,S.Steel SS304 2mm,RAL9016
- **HOLDER WALL-1200M-01**
Holder,Wall,1200M/S,S.Steel SS304 2mm,RAL9016

[1200 S/M](#)

[NR](#)



ADDED

• **NRP-USB-DOCKING-02**

NR Patch,USB Docking Station,NV,GmbH

• **BP-TANGO-M-STAND-01**

Deluxe Mobile Stand for tango requires Mount

• **BP-TANGO-MOUNT-01**

Mount for BP-TANGO-M-STAND-01

• **SE1200PRO-PAPERZFOLD-1**

Recording paper,Zfold

,210mmx140mmx120P,SE1200Pro

• **NECG301-PAPER-ROLL-01**

NECG-301,Thermal Recording paper,Roll,80mmx20m

• **NECG301-PAPER-ZFOLD-01**

NECG-301,Thermal Recording paper,Z-fold,80mmx70mmx200P



REMOVED

NSPIRO

- NSPIRO-01
- NSPIRO-02
- NSPIRO-AMBI-STICK-01
- NSPIRO-CALIB-PUMP-01
- NSPIRO-FK-01
- NSPIRO-FLOW-NOSE-01
- NSPIRO-SILICON-ADAP-01

PRICE UPDATES

• **CABLES**

C10 / C7 / C5 / C4 / C3



[Download the Distributors price list 2025 from the Distributors support portal](#)

