

Norav Medical Newsletter - February 2026



NEWSLETTER

As we conclude 2025, Norav Medical reflects on a year defined by meaningful progress, regulatory achievements, and continued innovation across our diagnostic portfolio. Throughout the year, we strengthened our compliance foundations, expanded our solutions in ECG and beyond, and introduced technologies designed to improve clinical workflows and diagnostic confidence. These milestones, achieved together with our global distributor network, position us strongly as we enter 2026, ready to further advance cardiac and pulmonary diagnostics worldwide.

*Product is important, but
having great product
distribution is more
important.*

LET'S MEET AT WHX DUBAI

At the upcoming WHX Dubai show, formerly Arab Health, Norav Medical will present its solutions from a centrally located 21-square-meter booth at a new, significantly larger exhibition venue. This year's expanded format offers an excellent opportunity to meet and explore our latest innovations. You can schedule



SCHEDULE A MEETING



Florida Pavilion, Booth #S11.F51

February 9–12, 2026 | Dubai Expo City, UAE



Medica 2025

Left to right: Arzit S.- Marcom Manager. Matthias K.– Sales Director Germany, Austria & Switzerland, Nario A.– VP Sales Middle East, Africa & CIS, Y.K. – Founder & CEO, Boris W. – Hardware R&D Engineering Manager, ordan M. – VP International Sales, Benny W.- Operations, Alex K.– Product & Customer Support Manager, JulioO.–Support, Gioele M.– GermanyOffice Manager,

Medica 2025

November 11-14, Düsseldorf

We had a highly successful presence at Medica 2025 in Düsseldorf, one of the world's leading medical trade fairs. We were proud to be among 4,821 exhibitors from 70 countries, with the event attracting 78,000 visitors from 160 nations, offering an exceptional platform for global connections and partnerships.

Heart Rhythm 2026

April 23 - 26 2026, Chicago

WE ARE EXCITED TO ANNOUNCE OUR FIRST-TIME PARTICIPATION AT HEART RHYTHM 2026 IN CHICAGO, ONE OF THE PREMIER GLOBAL EVENTS FOR CARDIAC ELECTROPHYSIOLOGY AND HEART

valued distributors and long-time customers while creating exciting new opportunities with potential partners in emerging markets. The face-to-face interactions and productive discussions reinforced the strength of our existing relationships and opened doors to new markets and possibilities.

WE LOOK FORWARD TO SEEING YOU AGAIN AT MEDICA 2026 IN DÜSSELDORF!

OUR PRESENCE IN THE NORTH AMERICAN MARKET, CONNECTING WITH LEADING CARDIOLOGISTS, ELECTROPHYSIOLOGISTS, AND HEALTHCARE PROFESSIONALS.

AT OUR BOOTH, WE WILL SHOWCASE OUR ADVANCED HOLTER AND REMOTE CARDIAC MONITORING SOLUTIONS, INCLUDING THE NR-314-P COMPACT 14-DAY HOLTER RECORDER AND TELEEVENT, OUR INNOVATIVE REMOTE ECG MONITORING SYSTEM DESIGNED FOR TELEMEDICINE AND HYBRID CARE MODELS. VISITORS WILL HAVE THE OPPORTUNITY TO SEE LIVE DEMONSTRATIONS, DISCUSS CLINICAL WORKFLOWS AND INTEGRATION SCENARIOS, AND EXPLORE DISTRIBUTION AND COLLABORATION OPPORTUNITIES WITH OUR TEAM.

VISIT US AT BOOTH #2137



ENHANCED 1200W



PRECISION YOU CAN TRUST. PERFORMANCE YOU CAN SELL.

The Enhanced 1200W is Norav Medical's flagship wireless stress ECG system, refined to meet the real-world demands of today's stress labs, hospitals, and sports cardiology centers.

Built on decades of clinical experience and trusted worldwide, the Enhanced 1200W delivers exceptional signal quality, reliable wireless performance, and improved usability - all in a compact, lightweight design optimized for daily clinical workflows.

READY TO BRING THE ENHANCED 1200W TO YOUR MARKET?

- ◆ REQUEST A QUOTE
- ◆ DISCUSS AVAILABILITY AND PRICING
- ◆ PLAN DEMO AND EARLY ORDERS

WHAT'S NEW - AND WHY IT MATTERS

• Lighter, More Comfortable Design

252 g housing with a slimmer profile

Improved patient comfort and easier handling during stress testing

• Upgraded Wireless Performance

New RF engine with +8 dB output

Extended wireless range and stronger connection stability, even under motion

• Enhanced Pacemaker Detection

High-speed digital pacemaker spike detection at 128,000 samples/second

Faster, more accurate recognition for confident interpretation in complex cases

• Cleaner ECG Under Motion

Refined signal processing and adaptive filtering

Clear, diagnostic-grade ECG traces during treadmill, stress echo, and CPET tests

• Simplified Power Management

Single AA battery design

Reduced weight, easier maintenance, and efficient operation

[LEARN MORE →](#)



NEW ACCESSORY: WALL HOLDERS

COMPATIBLE WITH ALL NR-SERIES AND 1200M/S DEVICES. A VALUABLE ADD-ON FOR YOUR CUSTOMERS SEEKING ORGANIZED, SPACE-EFFICIENT EQUIPMENT STORAGE. CONTACT US FOR DISTRIBUTOR PRICING AND AVAILABILITY.

STAY CONNECTED

FOLLOW US *in*

TANGO MOBILE

WIRELESS MOTION-TOLERANT BP & ECG TELEMETRY SYSTEM

COMING
SOON



COMING SOON

Get ready for Norav's breakthrough wearable monitoring system – a complete solution for real-time BP, ECG, and SpO₂ during patient movement. Tango Mobile is engineered to deliver accurate physiological data where traditional monitors fail – during walking, cycling, treadmill protocols, rehabilitation sessions, and mobile clinical care. This innovative system combines motion-tolerant blood pressure measurement, continuous 12-lead wireless ECG, and integrated pulse oximetry into one synchronized, body-worn platform. Designed for stress testing labs, cardiac and pulmonary rehabilitation centers, and sports medicine facilities,

THE NEW TANGO MOBILE IS COMING SOON – FIRST LIVE PREVIEW AT WHX DUBAI 2026. VISIT OUR BOOTH AND BE AMONG THE FIRST TO EXPERIENCE IT.

KNOW MORE:

● Motion-Tolerant Blood Pressure

ECG-gated auscultatory measurement ensures accurate BP readings during movement, with oscillometric mode for resting conditions.

● Wireless ECG Motion-Tolerant Telemetry

High-resolution, 12-lead wireless ECG with advanced motion-artifact reduction, delivering diagnostic-quality waveforms during physical activity.

● Integrated SpO₂ Monitoring

Continuous pulse oximetry adds essential respiratory data for complete cardiopulmonary assessment during exercise.

● Wearable Patient-Centered Design

Compact, lightweight, secure form factor enables natural movement without limiting patient activity or compromising signal quality.

● Unified Real-Time Display

Synchronized BP, HR, and SpO₂ values displayed on one integrated screen for immediate insight.

● Seamless System Integration

Wireless connectivity to NM-700 Central Station, stress ECG systems, and clinical workflows for comprehensive data management.

[LEARN MORE →](#)

SOFTWARE UPDATES

SOFTWARE
UPDATES


NEMS-Q (V2.7.9)

- **MiniSpir Integration** – full support for MiniSpir spirometer devices.
- **Manual Worklist Creation** – create worklists directly in NEMS-Q for 1200PRO device.

- **NR-314-P Device Preparation** – ECG verification for NR-314-P devices using GDT protocol.
- **Enhanced GDT Output** – improved GDT format for Holter report generation.
- **Updated HL7 Messaging** – latest HL7 message format standards implemented.


STRESS ECG (V5.988)

- **Separate Target HR Formulas** – independent heart rate target calculations for treadmill and ergometer protocols.
- **PCECG MSI Installation** – streamlined installation package for PC-ECG systems.
- **Vital Signs Settings Panel** – customizable color coding for all values, enlarged font for time, load, and protocol display.
- **Average Window Options** – toggle legend display off to maximize ST segment value visibility.



HOLTER NH-301 (V.4.1.0.0)

- **Extreme Events Management** – navigate, assign, and report extreme events (Fastest, Slowest, Longest, Shortest) with automatic flagging during analysis.
- **Enhanced Custom Reports** – additional fields and metrics for extreme events automatically included in summary templates.



EXTREME EVENTS MANAGEMENT

- **Automatic Extreme Event Flags** – system automatically identifies and flags the fastest, longest, slowest, and shortest events during analysis.
- **Manual Event Assignment** – right-click menu enables users to manually assign extreme event flags to specific events during review.
- **Visual Event Indicators** – extreme events are marked with dedicated icons in the top-right corner of event cells for instant recognition.
- **Quick Navigation Tool** – new lookup icon on the events toolbar enables instant jump to any flagged extreme event.

ENHANCED REPORTING & ANALYSIS

- **New Report Fields** – comprehensive extreme event data automatically included in generated reports.
- **Extreme Event Variables** – detailed metrics for lowest, highest, longest, and shortest cardiac events.
- **Mean Overall HR Metrics** – extended heart rate variables specifically calculated for extensive event analysis.

SYSTEM IMPROVEMENTS

- **Templates Folder Validation** – automatic verification of templates directory on first launch after installation, ensuring report generation integrity.
- **Single Flag Management** – intelligent flag system allows only one event per category, automatically transferring flags when reassigning.



NORAV CATALOG 2026

EXPLORE AND DOWNLOAD NORAV'S DETAILED AND COMPREHENSIVE PRODUCT CATALOG FROM THE DISTRIBUTOR PORTAL TODAY! IT'S AVAILABLE NOW AND READY TO SUPPORT YOUR BUSINESS.

[DOWNLOAD](#) →



2026 PRICING UPDATE

ORDER BY MARCH 31 TO
SECURE CURRENT PRICING

2026 PRICING UPDATE – APRIL 1 EFFECTIVE DATE

OUR 2026 DISTRIBUTOR PRICE LIST REFLECTS TARGETED ADJUSTMENTS ON TM6 ACCESSORIES, ERGOSELECT CYCLES, NBP-ONE ABPM SYSTEMS, AND TANGOM2 BP DEVICES, PLUS A PORTFOLIO-WIDE ADJUSTMENT TO ALIGN WITH CURRENT MARKET CONDITIONS.

ORDER BY MARCH 31, 2026 TO SECURE CURRENT PRICING.

REVIEW YOUR INVENTORY REQUIREMENTS NOW AND SUBMIT ORDERS TO MAXIMIZE PROFITABILITY BEFORE THE APRIL 1 TRANSITION. THE UPDATED PRICE LIST IS AVAILABLE IN OUR DISTRIBUTOR PORTAL.

NEW

Enhanced 1200W

CONTINUOUSLY
AVAILABLE

1200W





~~1200W E-[C10/S10]-19~~
1200W2-K-00-U-[C10/S10]-19

Wireless Stress,Kit

~~3. 1200W3-K-00-E-S10-02~~
1. 1200W3-K-00-E-C10-02

Wireless
Stress,Kit,NV,Clip,IEC,10L,WRV3,GmbH

~~4. 1200W3-K-00-U-S10-02~~
2. 1200W3-K-00-U-C10-02

Wireless
Stress,Kit,NV,Clip,AHA,10L,WRV3,GmbH

~~3. 1200W3-K-00-E-S10-02~~ ~~4. 1200W3-K-00-U-S10-02~~
Wireless Wireless
The Enhanced 1200W is now available alongside our standard 1200W model. Both versions remain in our active product line and are included in the updated price list.

Tango Mobile

1. With Nonin SPO2

SKU: **TANGO-MOB-ECG-KIT-12-[2/4]**

2. With Suntech SPO2

SKU: **TANGO-MOB-ECG-KIT-12-[1/3]**



DOWNLOAD →



info@norav.com | WWW.NORAV.COM

© 2026 Norav Medical Ltd. All rights reserved.
You are receiving this newsletter as a valued distributor or partner of Norav Medical.

[Unsubscribe](#) | [Update Preferences](#)



Follow US on LinkedIn



Copyright © 2026 Norav Medical, All rights reserved.

Want to change how you receive these emails?
You can [update your preferences](#) or [unsubscribe from this list](#).